

In the Claims:

Please cancel claims 2, 3, and 4.

Please amend claim 1 as follows.

1. (Currently Amended) A method ~~of using a set of calibration standards comprised of a plurality of ferromagnetic slugs~~ to provide a temperature calibration for a VTGA vacuum thermogravimetric analyzer (VTGA) comprising:
 - selecting a plurality of ferromagnetic slugs wherein each slug includes an alloy including nickel and copper and wherein the amount of copper is within the range of 15% to 28%;
 - determining the Curie temperature of a each slug;
 - placing the each slug in a sample holder of a VTGA within a magnetic field;
 - setting a temperature of the VTGA to a temperature corresponding to a first set-point temperature greater than the Curie temperature of the each slug by an amount equal to a first offset value;
 - holding the temperature of the VTGA at a temperature corresponding to a said first set-point temperature for a first time interval sufficient to allow the VTGA to thermally equilibrate;

increasing the temperature of the VTGA to a second set-point temperature greater than an immediately preceding set-point by an amount equal to a second offset value;

holding the temperature of the VTGA at a temperature corresponding to a said set-point temperature for a second time interval sufficient to allow the VTGA to thermally equilibrate; and,

if ~~the~~ each slug does not lose magnetization, then repeating the previous two operations, afterwards returning to this test; but if not, then recording the set-point temperature at which ~~the~~ each slug loses magnetization as the apparent Curie temperature of the each slug.

2. (Cancelled)

3. (Cancelled)

4. (Cancelled)